INSTRUCTION MANUAL

A Product by
Jyoti Ceramic Industries Pvt. Ltd.,(India)
Bio-Ceramic Division

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ABOUT JYODENT®

JYODENT® Blanks are manufactured from a clinical grade Yttria stabilized tetragonal zirconia polycrystal (Y-TZP) Zirconium Dioxide powder, pre-sintered to chalk hardness, suitable for manufacturing a wide variety of dental restoratives viz. coping, bridges, implant abutments etc that possess optimum bio compatibility, technical and aesthetic features.

JYODENT® Blanks are most gentle, provide chip free smooth milling of restoratives on open (not bar coded) CNC CAD-CAM milling Machines and MOM (Manually Operated milling machines) used by dental technicians.

Trademark	JYODENT®
Description	Oxide Ceramic Preparation in the form of blank for the manufacturing of dental frameworks
Material	Y-TZP Yttria (Y ₂ O ₃) stabilized tetragonal Zirconia (ZrO ₂) Polycrystal
Composition (in % by mass)	Zirconium Dioxide 94-96%, Yttrium Oxide 4-6%, Hafnium oxide 1-2%, Aluminium Oxide (Al ₂ O ₃) 0-0.5%

PHYSICAL & CHEMICAL PROPERTIES

PHYSICAL PROPERTIES

Colour	White	Ignition Point	Non Flammable
Surface Finish	Smooth	Pre-sintered Density	3.00 ± 0.2 g/cm ³
Condition	Solid Block	Final Sintered Density	6.08 ± 0.2 g/cm ³
Odour	None	Solubility after final sintering	Non Soluble
Melting Point	2680°C	Radioactivity	<0.01 Bq/g
Boiling Point	5500°C	(It has passed radioactivity test)	

CHEMICAL PROPERTIES

Ingredients	Percentage	CAS No.	EINECS No.
Zirconium Dioxide	94-96%	1314-23-4	215-227-2
Yttrium Oxide	4-6%	1314-36-9	215-233-5
Hafnium Oxide	1-2%	12055-23-1	235-013-2
Aluminum Oxide (Al ₂ O ₃)	0-0.5%	1344-28-1	215-691-6

Description : Porous ZrO₂ Block.

CERTIFICATIONS

JYODENT Zirconia Ceramic Dental Blanks are in compliance with the following standards

ISO 13356	Implants for surgery ceramic material based on Ytria stabilized Tetragonal Zirconia (Y-TZP)
ISO 13485	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 6872	Medical devices Dentistry in Ceramic Material.
ISO 7405	Dentistry Evaluation of Bio compatibility of Medical Devices used in dentistry.
MDR 2017	CDSCO Reg. No.: MFG/MD/2023/000011

HAZARDS AND THEIR EFFECTS

HAZARDOUS NATURE OF INGREDIENTS

(T : Total particulate matter)

Ingredients	ACGIH TLV	OSHA PEL	Units
Zirconia	5(T)	5(T)	mg/m ³
Yttria	2(T)	N.A.	mg/m ³

Hazard Identification

- Hazard Description : Not applicable
- Information concerning any particular hazards to man and environment : Not applicable.

Hazardous Nature of Ingredients

- The chemical ingredients are not hazardous to health. However, grinding or machining of this product causes dust and this dust may cause some irritation.
- Jyodent Zirconia Ceramic Dental Blanks do not contain any hazardous substances, specifically phthalates.
- There is no leakage from the device.

First AID measures

- Inhalation** : In case dust inhaled and there are symptoms of cough or shortness of breath, move affected person away from exposure area and into a well ventilated place. If breathlessness recurs, seek medical attention.
- Skin Contact** : If irritation to skin occurs due to contact with dust then wash the affected area with clean fresh water. If irritation persists, seek medical attention.
- Eye Contact** : If the dust comes in contact with the eyes and irritation to eyes occurs, wash the affected eye with fresh clean tap water.
- Ingestion** : If substantial quantity of dust is ingested, drink large quantity of water and seek medical attention.

Fire Safety Measures:

- No special fire safety measures are required.
- This product is fire resistant, non-combustible and non-hazardous.

PROTECTIVE MEASURES

Accidental Release Measures

- Personal Safety Precautions : Use protective clothing, eye gear and dust mask while machining of blocks.

Environmental Protection

- No special measures required.

Cleaning Methods

- Collect the dust and pieces by vacuum or mechanically.
- Land filling to be done as per regulations set by the local authorities.

Handling

- As per Instruction Manual. Use protective clothing, eye gear and dust mask while machining of blanks.
- No other special measures necessary.

Storage

- Material should be stored in well-protected, cool and dry place in its original packaging.

Exposure

- Avoid exposure and maintain good housekeeping procedures.
- Ensure adequate ventilation and suitable dust extraction system in the work place.

Reactivity Data

- Material is not reactive under normal conditions.
- Material is stable with no hazardous decomposition.

Toxicological Information

- There are no known reactions which are damaging to health.

Ecological Information

- Do not mix/empty the product or collected dust into water or drains.

Disposal Considerations

- Product disposal must be done as per regulations specified by the local authorities.
- Packaging disposal must be done as per regulations specified by the local authorities.

HANDLING PROCESS

Handling:

JYODENT Blanks can be milled on any dental open(non-barcoded) CNC CAD-CAM machines as well as on manually operated copy milling machines (MOM) to procure Zirconia copings, Bridges & other full ceramic restoratives.

Direction for use:

- 1) Jyodent blanks should be milled to the required frame shape. The milling must be approximately 25% to the original size to compensate for shrinkage during the final sintering process.
- 2) The exact shrinkage percentage or expansion factor should be programmed into the Cam software. Kindly refer the exact individual shrinkage factor mentioned on each surface of the Jyodent discs & blanks as well as on Jyodent product manuals.

GRINDING & POLISHING

- 1) After milling the frame work, It can be further shaped with a hand piece & diamond bur. Gently, at under 10,000 rpm. The frame work must always be wet ground to avoid overheating.
- 2) It is recommended to sand blast with 50-micron aluminum oxide at 2-4 bar, it clean the framework & recover any strength loss during milling.

Alternatively


- 3) After the final sintering the zirconia restorative can be fitted to the model by micro motor* with diamond coated burrs. Over heating has to be avoided.
- 4) Finally clean the framework using ultrasonic cleaner or steam cleaner.

SINTERING PROCESS

1. Firstly, inspect the colour, shape, number of Duralox alumina containing (99.7% Al_2O_3) ceramic beads. If yellowish and mutilated ceramic beads are observed they should be replaced immediately.
2. Inspect the sintering furnace. It is important to clean the furnace cavity on time.
3. Lay the JYODENT restorative on loose high purity (99.7% Al_2O_3) ceramic micro beads loaded in a refractory ceramic saggar.
4. Apply gentle pressure on the restorative so that they take proper support on the loose ceramic beads.
5. Ensure that no ceramic bead gets jammed in the internal cavities of the framework.
6. For long span bridges, kindly keep the supporting framework on the lingual side to prevent distortion or rocking of Zirconia bridge post sintering.

SINTERING PROGRAMME FOR JYODENT RESTORATIVES

- For Coping and bridges : From ambient ($25^\circ\text{C}/77^\circ\text{F}$) to final sintering temperature ($1450^\circ\text{C}/2642^\circ\text{F}$) heating rate $8^\circ\text{C}/\text{min}$. Soak/Hold period in kiln 2.0 hrs.
- Let the kiln cool down to ambient temperature naturally before opening.

Note : Heating & Cooling rate depends upon furnace capacity.

SHRINKAGE CHART

Please refer following Shrinkage Chart before milling Jyodent Ceramic Dental Discs/Blanks.

Note : Shrinkage percentage is printed on surface of each Jyodent Ceramic Dental Disc & Blank.

Shrinkage %	Scaling factor	Shrinkage %	Scaling factor
18.00%	1.219	18.10%	1.221
18.20%	1.222	18.30%	1.223
18.40%	1.225	18.50%	1.226
18.60%	1.228	18.70%	1.230
18.80%	1.231	18.90%	1.233
19.00%	1.234	19.10%	1.236
19.20%	1.237	19.30%	1.239
19.40%	1.240	19.50%	1.242
19.60%	1.243	19.70%	1.245
19.80%	1.246	19.90%	1.248
20.00%	1.250	20.10%	1.251

Shrinkage %	Scaling factor
20.20%	1.253
20.40%	1.256
20.60%	1.259
20.80%	1.262
21.00%	1.265
21.20%	1.269
21.40%	1.272
21.60%	1.275
21.80%	1.278
22.00%	1.282

Shrinkage %	Scaling factor
20.30%	1.254
20.50%	1.257
20.70%	1.261
20.90%	1.264
21.10%	1.267
21.30%	1.270
21.50%	1.273
21.70%	1.277
21.90%	1.280

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Jyoti Ceramic GmbH
Techno Ceramic Inc.